

Evaluation of a Nitrogen Flush System to Prevent Oxidation of Fish Oil Encapsulated in Licaps® Capsules Using CFS1200 Equipment

S. Robin¹, S. Grunenwald², P. Peter³, J-L. Colin², D. Cade¹

¹Chemical R&D Department, Capsugel, Division of Pfizer,

² Pharmaceutical R&D Center, Capsugel, Division of Pfizer,

³ Pharmaceutical Technology Group, Capsugel, Division of Pfizer


¹Dominique.Cade@pfizer.com

Poster presented at the 2009 Annual Meeting and Exposition of the American Association of Pharmaceutical Scientists.

Los Angeles, California

November 8-12, 2009

BAS 411



Evaluation of a Nitrogen Flush System to Prevent Oxidation of Fish Oil Encapsulated in Licaps® Capsules Using CFS1200 Equipment

S. Robin ¹, S. Grunenwald ², P. Peter ³, J-L. Colin ², D. Cade ¹

¹Chemical R&D Department, Capsugel, Division of Pfizer,

² Pharmaceutical R&D Center, Capsugel, Division of Pfizer,

³ Pharmaceutical Technology Group, Capsugel, Division of Pfizer

Key words: encapsulation, oxidation prevention, nitrogen, capsules, Fish oil, stability

PURPOSE

The work presented hereafter demonstrates the efficiency of the Nitrogen Flush system installed on CFS1200 machine (Photo 1) to reduce oxidation of oil encapsulated in Licaps® capsules.




Photo 2: Process of encapsulation on CFS1200 machine



Photo 1: CFS1200 machine


Unsaturated oils are prompt to oxidation in contact with oxygen from the air (Fig. 1). During encapsulation in hard capsules, these oils are in contact with the air at filling and after the closing steps (Photo 2, orange arrows at filling and after closing).

This air can be replaced by nitrogen or another inert gas before closing thanks to Nitrogen Flush accessories. The stability of the oil encapsulated in these conditions increases to a large extent. The purpose of this study is to evaluate how long the capsules will remain below the limit of oxidation when stored under I.C.H. conditions.

METHODS

Product description

The capsules considered are Licaps® capsules size #0, natural transparent (Capsugel). The fish oil encapsulated is EPAX6000 TG (Polaris). The nitrogen used is technical grade with less than 5ppm traces of oxygen (Gaz Liquide). The sealing fluid is a mixture of pure ethanol and water 50/50 w/w.



Encapsulation conditions

- A design of experiments was performed using two factors at two levels: nitrogen flush (Yes/No) and product degassing in the CFS1200 container prior to encapsulation (Yes/No). The four runs were repeated the day after to confirm repeatability of results.
- Each Licaps® capsule was filled with 500mg fish oil, closed and sealed in the presence of 20µL sealing fluid on CFS1200. Nitrogen flow was 11L/min at closing, and 2L/min during 1 hour for degassing the oil in the CFS1200 container.
- Samples of oil were taken before and after encapsulation of each run; samples of filled capsules were stored in HDPE bottles in climatic cabinets under I.C.H. conditions (25°C/60%RH, 30°C/65%RH and 40°C/75%RH) for 6 months to study their stability over time in standard and accelerated conditions.

Measurement of oil oxidation

The oxidation of oil is measured using two parameters:

- Peroxide index (IP): milliequivalents of peroxide per kg of fat are measured by titration with iodide ion.
- Anisidine value (AV): this method measures the level of non-volatile carbonyl components (aldehydes and ketones), which are formed during deterioration of oils. They react with p-anisidine determining an absorption that can be measured at 350 nm.

According to European Pharmacopeia, maximum limits for oxidation products are IP<10 and AV<30. IP points out the first phase of oxidation, whereas AV stands for the second phase of oxidation that ends up with total degradation of oil (brownish color and release of rancid odor) (Fig 1).





Fig. 1: Oxidation of unsaturated oils: (A) unsaturations, (B) peroxides, (C) products of degradation, (D) volatile compounds (odor)

RESULTS

The overall results are consistent between repeats (Fig 2: day factor).

Nitrogen flush is the main significant factor that affects both IP and AV: the oxidation values are lower with nitrogen flush.

The fact of degassing the oil in the container prior to encapsulation is not significant in our processing conditions (4 hours running per day).






Fig. 2: Main effects of factors “nitrogen flush”, “degassing” and “day”

The detailed results show that peroxide index (IP) increases significantly after encapsulation when nitrogen flush is not used (Fig. 3). The IP value is larger in standard conditions (25°C/60%RH) than in accelerated ones (40°C/75%RH) because peroxides destroy each other beyond a threshold of concentration, as described in literature.

Anisidine value (AV) consistently increases over the 6 months whatever the conditions. As expected, the oxidation is much more rapid under accelerated conditions: AV limit is reached after 6 months stability in oil encapsulated without nitrogen flush.




Fig. 3: Evolution of peroxide index oil encapsulated over time

The results are conforming to the ones expected (Table 1): fish oil oxidation is significantly slowed down thanks to nitrogen flush during closing.

Under accelerated conditions, Fish oil oxidation is prevented for 6 months with nitrogen flush: this corresponds to 2 years shelf-life extension under standard conditions.


Degassing of oil in CFS1200 container has minor effect on oxidation because encapsulation time was short; it would have been significant at larger scale.

Runs	Nitrogen flush	Oil degassing	T=0	1 Month	3 Months	6 Months
4 & 8	No	No	IP>10	IP>10	IP>10	AV>30 @40°C/75%RH
2 & 6	No	Yes	IP>10	IP>10	IP>10	AV>30 @40°C/75%RH
3 & 7	Yes	No	OK	OK 3C	OK 3C	OK 3C
1 & 5	Yes	Yes	OK	OK 3C	OK 3C	OK 3C

Table 1: Stability of oil encapsulated. “OK 3C”:IP<10 and AV<30 at the 3 ICH conditions

CONCLUSIONS

Fish oil encapsulated in Licaps® capsules under nitrogen is stable after 6 months storage under standard and accelerated conditions thanks to nitrogen flush. This means the product has 2 years shelf life instead of 1 year without nitrogen flush.



REFERENCES

- Przybylski R. et al., J. Am. Oil Chem. Soc. (1988), 65:629
- Labuza T.P., CRC Crit. Rev. Food Technol. (1971), 2:355
- Chan, H.W.S “autoxidation of unsaturated Lipids” (1987, Ed. H.W.S.Chan, Academic press, New York;
- J.K. Kaitaranta, J. Am. Oil Chem. Soc., (1992), Vol. 69, no. 8
- P.J.Ke at al., Spring Am. Oil Chem. Soc. Meeting, (1974) Mexico City.
- CFS1200™ brochure references
- I.C.H. Q1A(R2): Stability testing of new drug substances and products
- International Stability Testing (1999), D.J. Mazzo, Interpharm Press

ACKNOWLEDGEMENT

We wish to acknowledge the teams from the Chemical R&D, the Analytical R&D, the Pharmaceutical Technology Group and Dr. Keith Hutchison for their support and discussions.

Poster presented at the 2009 Annual Meeting and Exposition of the American Association of Pharmaceutical Scientists.
Los Angeles, California
November 8-12, 2009

201006002

