



Capsule Count Guide

Capsugel offers various packing options to ensure you get the quantity of capsules required in the best condition. Storage and handling information and recommendations are located on the outside of all packaging options as well as in the Capsugel Technical Reference Files. To ensure optimum capsule performance, please follow these recommendations.

Cartons

In order to protect Capsugel capsules from large variations in the surrounding relative humidity, standard packaging consists of protective, food grade antistatic bags inside robust cardboard cartons. Each bag is tamper-evident sealed ensuring integrity of the product. Capsugel also offers a fiber-free carton option made of type 5 plastic which can be recycled to protect capsules.

Dimensions

Cardboard Case & Fiber-Free Case 23-5/8" L x 15-5/8" W x 28-3/4" H

Capsule Caddy™ Container

Made from 100% recyclable, high-density polyethylene to protect capsules from the degradations of light and moisture.

Dimensions

Capsule Caddy[™] 11-15/16" L x 9-7/8" W x 14-11/16" H

Small Box

12" L \times 9-15/16" W \times 9-3/4" H

R&D Jars

The R&D jars are made from white HDPE plastic with a lid. HDPE is a Type 2 plastic suitable for recycling.

Dimensions

Height: 8-5/16" Diameter: 4-3/4"

Notes

- 2. The Capsule Caddy, R&D Jar, and Small Box have +/- 10% quantity tolerance.



Capsule Count Guide

| SIZE | QUANTITY PER CARTON | QUANTITY PER CAPSULE CADDY | QUANTITY PER R&D JAR | | |
|--|----------------------------------|-------------------------------|-------------------------|--|--|
| Coni-Snap® | Coni-Snap® Hard Gelatin Capsules | | | | |
| 000 | 50,000 | 6,500 | 800 | | |
| 00E | 65,000 | 9,000 | 1,075 | | |
| 00 | 70,000 | 9,000 | 1,250 | | |
| OEL | 90,000 | 11,000 | 1,500 | | |
| 0 | 100,000 | 12,000 | 1,500 | | |
| 1E | 115,000 | 16,500 | 1,900 | | |
| 1 | 125,000 | 17,000 | 2,300 | | |
| 2EL | 150,000 | 15,000 | 2,500 | | |
| 2 | 175,000 | 21,000 | 3,200 | | |
| 3 | 225,000 | 28,500 | 3,800 | | |
| 4 | 300,000 | 36,500 | 5,000 | | |
| 5 | 350,000 | 50,000 | 5,000 | | |
| Coni-Snap® | Sprinkle Capsules | - 5 | | | |
| 0* | 100,000 | 12,000 | 1,500 | | |
| Licaps® Capsules (indicate banding or LEMS® Sealing) | | | | | |
| 00* | 70,000 | 9,000 | 1,250 | | |
| 0* | 100,000 | 12,000 | 1,500 | | |
| 1 | 125,000 | 17,000 | 2,300 | | |
| 2 | 175,000 | 21,000 | 3,200 | | |
| 3 | 225,000 | 28,500 | 3,800 | | |
| Vcaps® Capsules (hypromellose) | | | | | |
| 000 | 50,000 | 6,500 | 800 | | |
| 00E | 65,000 | 9,000 | 1,075 | | |
| 00 | 70,000 | 9,000 | 1,075 | | |
| OEL | 90,000 | 11,000 | 1,500 | | |
| 0 | 100,000 | 12,000 | 1,500 | | |
| 1 | 125,000 | 17,000 | 2,300 | | |
| 2 | 175,000 | 21,000 | 3,200 | | |
| 3 | 225,000 | 28,500 | 3,800 | | |
| 4 | 300,000 | 36,500 | 3,800 | | |
| Vcaps® Plus Capsules (hypromellose) | | | | | |
| 00E | 65,000 | 9,000 | 1,075 | | |
| 00 | 70,000 | 9,000 | 1,250 | | |
| OEL | 90,000 | 11,000 | 1,500 | | |
| 0 | 00,000 | 12,000 | 1,500 | | |
| 1E | 115,000 | N/A | N/A | | |
| 1 | 125,000 | 17,000 | 2,300 | | |
| 2 | 175,000 | 21,000 | 3,200 | | |
| 3 | 225,000 | 28,500 | 3,800 | | |
| 4 | 300,000 | 36,500 | 5,000 | | |
| 5 | 350,000 | 50,000 | 6,500 | | |

| SIZE | QUANTITY PER CARTON | QUANTITY PER CAPSULE CADDY | QUANTITY PER R&D JAR | | |
|-------------------------------------|--------------------------------|--|-------------------------|--|--|
| All Color Capsules | | | | | |
| 0 | 100,000 | 12,000 | N/A | | |
| 1 | 125,000 | 17,000 | N/A | | |
| 2 | 175,000 | 21,000 | N/A | | |
| 3 | 225,000 | 28,500 | N/A | | |
| 4 | 300,000 | 36,500 | N/A | | |
| DRcaps™ Acid-Resistant Capsules | | | | | |
| 00 | 70,000 | Please contact your Customer Service Representative for | 1,250 | | |
| 0 | 100,000 | | 1,500 | | |
| 1 | 125,000 | | 2,300 | | |
| 2 | 175,000 | | 3,200 | | |
| 3 | 225,000 | availability | N/A | | |
| Plantcaps® | Plantcaps® Capsules (pullulan) | | | | |
| 00 | 70,000 | Please contact your Customer Service Representative for availability | 1,250 | | |
| 0 | 100,000 | | 1,500 | | |
| 1 | 125,000 | | 2,300 | | |
| 2 | 175,000 | | N/A | | |
| 3 | 225,000 | | N/A | | |
| 4 | 300,000 | | N/A | | |
| DPI®, Vcap® or Vcap Plus® Capsules | | | | | |
| 3 | 225,000 | 28,500 | 3,800 | | |
| DBcaps® Over-encapsulation Capsules | | | | | |
| AAA | 50,000 | 5,500 | 2,500 | | |
| AAE* | 50,000 | 7,000 | 3,800 | | |
| AA* | 50,000 | 9,500 | 4,700 | | |
| A* | 90,000 | 12,000 | 5,000 | | |
| В | 125,000 | 17,000 | 8,000 | | |
| С | 150,000 | 21,800 | 9,500 | | |
| D | 200,000 | 33,000 | 11,000 | | |
| Е | 250,000 | N/A | N/A | | |

^{*} Available in gelatin and Vcaps Plus.

For more information

Visit Capsugel.com or contact your Capsugel representative at 800-845-6973.